



More information on the website  
[mirror.radwag.com/it/info,w1,D5Z](https://mirror.radwag.com/it/info,w1,D5Z)

# Cilindrico, set (1 mg - 200 g), cassa di legno

OK-501-0039



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

Parametri metrologici	
Classe OIML	E1
Costruzione	
Materiale	Stainless steel
Peso	611,11 g
Valore nominale	1 mg – 200 g
Numero	23 pcs

## Compatibile con (Additional Fee)

Pinzette

Pennelli