



More information on the website
mirror.radwag.com/de/info,w1,D64

Prüfgewichte der M1 Klasse - Knopfgewichte, Set (1 mg - 200 g), Kunststoff-Etui

OK-402-0073



The drawings, photos and graphics used are for illustrative purposes only.

Technische Daten

Messtechnische Parameter	
OIML-Klasse	M1
Gewicht	611,11 g
Nominalwert	1 mg – 200 g
Anzahl	23 St.

Kompatibel mit (Additional Fee)

Pinsel