



More information on the website  
[mirror.radwag.com/en/info,w1,DBV](https://mirror.radwag.com/en/info,w1,DBV)

# M1 Mass Standard - Knob Weights - 200 kg

OP-403-0109



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

### Metrological parameters

OIML Class	M1
------------	----

### Construction

Material	Stainless steel
----------	-----------------

Nominal value	200 kg
---------------	--------

## Compatible with (Additional Fee)

Brushes