



More information on the website
mirror.radwag.com/en/info,w1,DQ9

BT 300 μ l Tips - 1000 pcs

PT-111-0001



The drawings, photos and graphics used are for illustrative purposes only.

Overview

OMNITIP™ BT 300 μ l tips for pipettes are high quality laboratory accessories designed for precision pipetting. Designed with efficiency and convenience in mind, they are an indispensable tool for everyday laboratory work.

Capacity: 300 μ l Colour: transparent Set contains 1000 tips in a plastic bag with a zipper closure

Application: - precise dispensing of biological samples, - accurate pipetting of medicines and medical solutions, - preparing and transferring chemical solutions, - an indispensable tool in experiments requiring precise dosing.

Compatible with (Additional Fee)

Automatic Variable-Volume Pipettes

Workstation for Pipettes Calibration