



More information on the website
mirror.radwag.com/en/info,w1,DQ9

BT 300 µl Tips - 1000 pcs

PT-111-0001



The drawings, photos and graphics used are for illustrative purposes only.

Overview

OMNITIP™ BT 300 µl tips for pipettes are high quality laboratory accessories designed for precision pipetting. Designed with efficiency and convenience in mind, they are an indispensable tool for everyday laboratory work.

Capacity: 300 µl Colour: transparent Set contains 1000 tips in a plastic bag with a zipper closure

Application: - precise dispensing of biological samples, - accurate pipetting of medicines and medical solutions, - preparing and transferring chemical solutions, - an indispensable tool in experiments requiring precise dosing.

Compatible with (Additional Fee)

Automatic Variable-Volume Pipettes

Workstation for Pipettes Calibration