



More information on the website
mirror.radwag.com/en/info,w1,DQP

F2 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 mg - 50 g), Plastic Box

OK-302-0034



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Metrological parameters | |
|-------------------------|-----------------|
| OIML Class | F2 |
| Construction | |
| Material | Stainless steel |
| Weight | 111.11 g |
| Nominal value | 1 mg – 50 g |
| Quantity | 20 pcs |

Compatible with (Additional Fee)

Tweezers

Brushes