



More information on the website  
[mirror.radwag.com/en/info,w1,DVQ](http://mirror.radwag.com/en/info,w1,DVQ)

# F1 Mass Standard - Knob Weights Without Adjustment Chamber, Set (1 mg - 10 kg), Plastic Box

OK-202-0088



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F1
Construction	
Material	Stainless steel
Weight	21 111.11 g
Nominal value	1 mg – 10 kg
Quantity	29 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes