



More information on the website  
[mirror.radwag.com/it/info,w1,DW8](https://mirror.radwag.com/it/info,w1,DW8)

## Cilindrico, set (1 mg - 500 g), imballaggi in plastica



The drawings, photos and graphics used are for illustrative purposes only.

### Specifiche

Parametri metrologici	
Classe OIML	F1
Costruzione	
Materiale	Stainless steel
Peso	1 111,11 g
Valore nominale	1 mg - 500 g
Numero	24 pcs

### Compatibile con

Pinzette

Pennelli