



More information on the website
mirror.radwag.com/it/info,w1,DXV

Valigia di trasporto per Termobilancia MA X2, MAR

WX-003-0023



The drawings, photos and graphics used are for illustrative purposes only.

Specifiche

Parametri fisici

Dimensioni del pacco	545×330×340 mm
----------------------	----------------

Compatibile con (Additional Fee)

Termobilancia MA X7.IC.A
Termobilancia MA 5Y
Termobilancia MA X2.IC.A
Termobilancia MA X2.A

Termobilancia MA X7.A
Termobilancia MA 5Y.IC
Termobilancia MA R

Device dimensions

