



More information on the website
mirror.radwag.com/en/info,w1,E8Z

F2 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 g - 500 g), Wooden Box

OK-301-0030



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Metrological parameters | |
|-------------------------|-----------------|
| OIML Class | F2 |
| Construction | |
| Material | Stainless steel |
| Weight | 1 110 g |
| Nominal value | 1 – 500 g |
| Quantity | 12 pcs |

Compatible with (Additional Fee)

Tweezers

Brushes