



More information on the website
mirror.radwag.com/it/info,w1,EIB

Cilindrico, set (1 mg - 5 kg), imballaggi in plastica

OK-202-0087



The drawings, photos and graphics used are for illustrative purposes only.

Specifiche

| Parametri metrologici | |
|-----------------------|-----------------|
| Classe OIML | F1 |
| Costruzione | |
| Materiale | Stainless steel |
| Peso | 11 111,11 g |
| Valore nominale | 1 mg - 5 kg |
| Numero | 28 pcs |

Compatibile con (Additional Fee)

Pinzette

Pennelli