



More information on the website  
[mirror.radwag.com/en/info,w1,EVV](https://mirror.radwag.com/en/info,w1,EVV)

# Knob weights, Set (1 mg - 1 kg), Wooden Box



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	E2
Construction	
Material	Stainless steel
Weight	2 111.11 g
Nominal value	1 mg – 1 kg
Quantity	25 pcs

## Compatible with

Tweezers

Brushes