



More information on the website  
[mirror.radwag.com/br/info,w1,F64](https://mirror.radwag.com/br/info,w1,F64)

# Cable to the ZR-03/1 and ZR-04/1 Powerbank

WX-013-0244



The drawings, photos and graphics used are for illustrative purposes only.

## Visão geral

It is 20 cm long. It is compatible with: 5Y, 4Y, X2, C32, R2, WLC, WTC, C315, CY10, WLY and WPY series balances and scales.

## Compatível com (Additional Fee)

Microbalança para calibração de micropipetas XA 5Y.M.A.P  
 Plataforma multifuncional WPY  
 Microbalança para calibração de micropipetas MYA 5Y.P  
 Plataforma de pesagem C315  
 Balança de precisão 5Y.PM  
 Microbalança para filtros MYA 5Y.F  
 Balança analítica AS X7  
 Microbalança XA 5Y.M.A  
 Balança de precisão WLC  
 Balança de precisão PS R2  
 Balança de precisão WLC X7  
 Terminal de pesagem PUE CY10  
 Plataforma de pesagem C315.4

Balança analítica AS 5Y  
 Balança de precisão PS R2.H  
 Balança de precisão WTC  
 Balança de precisão CY10  
 Balança analítica para filtros XA 5Y.F  
 Balança de precisão WLY  
 Balança de precisão WLC C/2  
 Microbalança MYA 5Y  
 Balança de precisão PS R1  
 Balança analítica XA 5Y  
 Balança analítica AS R2 PLUS  
 Balança de precisão PS 5Y  
 Microbalança XA 5Y.M

Balança de precisão C32.PM  
Balança de precisão PS X7  
Balança analítica XA 5Y.A  
Ultra-Microbalança UYA 5Y

Balança analítica AS R1 PLUS  
Ultra-Microbalança para filtros UYA 5Y. F  
Microbalança para pesagem de stents XA 5Y.M.A.S  
Impressora de recibos