



More information on the website
mirror.radwag.com/br/info,w1,F64

Cable to the ZR-03/1 and ZR-04/1 Powerbank

WX-013-0244



The drawings, photos and graphics used are for illustrative purposes only.

Visão geral

It is 20 cm long. It is compatible with: 5Y, 4Y, X2, C32, R2, WLC, WTC, C315, CY10, WLY and WPY series balances and scales.

Compatível com (Additional Fee)

Microbalança para calibração de micropipetas XA 5Y.M.A.P
Plataforma multifuncional WPY
Microbalança para calibração de micropipetas MYA 5Y.P
Plataforma de pesagem C315
Balança de precisão 5Y.PM
Microbalança para filtros MYA 5Y.F
Balança analítica AS X7
Microbalança XA 5Y.M.A
Balança de precisão WLC

Balança de precisão PS 5Y
Balança de precisão PS R2.H
Balança de precisão WTC
Balança de precisão CY10
Balança analítica para filtros XA 5Y.F
Balança de precisão WLY
Balança de precisão PS X7
Balança de precisão WLC C/2
Microbalança MYA 5Y

Balança de precisão PS R2
Balança de precisão WLC X7
!Terminal de pesagem PUE CY10
Plataforma de pesagem C315.4
Balança de precisão C32.PM
Balança analítica XA 5Y.A
Ultra-Microbalança UYA 5Y
Balança analítica AS 5Y

Balança de precisão PS R1
Balança analítica XA 5Y
Balança analítica AS R2 PLUS
Microbalança XA 5Y.M
Balança analítica AS R1 PLUS
Ultra-Microbalança para filtros UYA 5Y. F
Microbalança para pesagem de stents XA 5Y.M.A.S
Impressora de recibos