



More information on the website
mirror.radwag.com/br/info,w1,F64

Cable to the ZR-03/1 and ZR-04/1 Powerbank

WX-013-0244



The drawings, photos and graphics used are for illustrative purposes only.

Visão geral

It is 20 cm long. It is compatible with: 5Y, 4Y, X2, C32, R2, WLC, WTC, C315, CY10, WLY and WPY series balances and scales.

Compatível com (Additional Fee)

Microbalança para calibração de micropipetas XA 5Y.M.A.P
 Plataforma multifuncional WPY
 Microbalança para calibração de micropipetas MYA 5Y.P
 Plataforma de pesagem C315
 Balança de precisão 5Y.PM
 Microbalança para filtros MYA 5Y.F
 Balança analítica AS X7
 Microbalança XA 5Y.M.A
 Balança de precisão WLC
 Balança de precisão PS R2
 Balança de precisão WLC X7
 Terminal de pesagem PUE CY10
 Plataforma de pesagem C315.4

Balança analítica AS 5Y
 Balança de precisão PS R2.H
 Balança de precisão WTC
 Balança de precisão CY10
 Balança analítica para filtros XA 5Y.F
 Balança de precisão WLY
 Balança de precisão WLC C/2
 Microbalança MYA 5Y
 Balança de precisão PS R1
 Balança analítica XA 5Y
 Balança analítica AS R2 PLUS
 Balança de precisão PS 5Y
 Microbalança XA 5Y.M

Balança de precisão C32.PM
Balança de precisão PS X7
Balança analítica XA 5Y.A
Ultra-Microbalança UYA 5Y

Balança analítica AS R1 PLUS
Ultra-Microbalança para filtros UYA 5Y. F
Microbalança para pesagem de stents XA 5Y.M.A.S
Impressora de recibos