



More information on the website  
[mirror.radwag.com/us/info,w1,F64](https://mirror.radwag.com/us/info,w1,F64)

## Cable to the ZR-03/1 and ZR-04/1 Powerbank



The drawings, photos and graphics used are for illustrative purposes only.

### Overview

It is 20 cm long. It is compatible with: 5Y, 4Y, X2, C32, R2, WLC, WTC, C315, CY10, WLY and WPY series balances and scales.

### Compatible with

XA 5Y.M.A.P Microbalance for pipettes calibration  
 WPY Multifunctional Scale  
 MYA 5Y.P Microbalance for pipettes calibration  
 Tunnel detectors  
 C315 Platform scale  
 5Y.PM Precision balance  
 MYA 5Y.F Microbalance for Filters  
 AS X7 Analytical Balance  
 PS X2 Precision Balance  
 AS X2 PLUS Analytical Balances  
 C315.D Baby scale  
 XA 5Y.M.A Microbalance  
 WLC Precision Balance  
 PS R2 Precision Balance

PS 5Y Balances for PGC  
 PS R2.H Precision Balance  
 WTC Precision Balance  
 CY10 Precision balance  
 C315.OW Personal scale  
 Analytical Balance for Filters XA 5Y.F  
 WLY Precision Balance  
 WLC C/2 Precision balance  
 MYA 5Y Microbalance  
 PS R1 Precision Balance  
 XA 5Y Analytical Balance  
 AS R2 PLUS Analytical Balance  
 PUE C32 Weighing Terminal  
 PS 5Y Precision Balance

WLC X7 Precision balance  
PUE CY10 Weighing Terminal  
C315.4 Platform Scale  
C32.PM Precision Balance  
C315.O Personal scale  
PS X7 Precision Balance  
XA 5Y.A Analytical Balance  
C315.4P pallet scales  
UYA 5Y Ultra-Microbalance  
AS 5Y Analytical Balance

XA 5Y.M Microbalance  
PUE C315 Terminal  
AS R1 PLUS Analytical Balance  
UYA 5Y.F Ultra-Microbalance for Filters  
C315.4B Bed scale  
C315.K Chair scale  
XA 5Y.M.A.S Microbalance for stents weighing  
C315.8B Ramps bed scale  
Receipt Printer