



More information on the website  
[mirror.radwag.com/en/info,w1,FEZ](https://mirror.radwag.com/en/info,w1,FEZ)

# F1 Mass Standard - Knob Weights With Adjustment Chamber - 20 g

OP-203-0088



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

### Metrological parameters

OIML Class	F1
------------	----

### Construction

Material	Stainless steel
----------	-----------------

Nominal value	20 g
---------------	------

## Compatible with (Additional Fee)

Tweezers

Brushes