



More information on the website  
[mirror.radwag.com/en/info,w1,F07](http://mirror.radwag.com/en/info,w1,F07)

# F2 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 g - 10 kg), Wooden Box

OK-301-0034



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F2
Construction	
Material	Stainless steel
Weight	21 110 g
Nominal value	1 g – 10 kg
Quantity	17 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes