



More information on the website  
[mirror.radwag.com/en/info,w1,FRE](http://mirror.radwag.com/en/info,w1,FRE)

# F1 Mass Standard - Knob Weights Without Adjustment Chamber, set (1 mg - 2 kg), Plastic Box

OK-202-0086



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | F1              |
| Construction            |                 |
| Material                | Stainless steel |
| Weight                  | 6 111.11 g      |
| Nominal value           | 1 mg – 2 kg     |
| Quantity                | 27 pcs          |

## Compatible with (Additional Fee)

Tweezers

Brushes