



More information on the website
mirror.radwag.com/en/info,w1,FRE

F1 Mass Standard - Knob Weights Without Adjustment Chamber, set (1 mg - 2 kg), Plastic Box



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Metrological parameters | |
|-------------------------|-----------------|
| OIML Class | F1 |
| Construction | |
| Material | Stainless steel |
| Weight | 6 111.11 g |
| Nominal value | 1 mg – 2 kg |
| Quantity | 27 pcs |

Compatible with

Tweezers

Brushes