



More information on the website
mirror.radwag.com/en/info,w1,G2P

F2 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 mg - 5 kg), Wooden Box

OK-301-0043



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Metrological parameters	
OIML Class	F2
Construction	
Material	Stainless steel
Weight	11 111.11 g
Nominal value	1 mg – 5 kg
Quantity	28 pcs

Compatible with (Additional Fee)

Tweezers

Brushes