



## Cilindrico, set (1 mg - 5 kg), cassa di legno

OK-301-0043

More information on the website  
[mirror.radwag.com/it/info,w1,G2P](http://mirror.radwag.com/it/info,w1,G2P)



The drawings, photos and graphics used are for illustrative purposes only.

### Specifiche

Parametri metrologici	
Classe OIML	F2
Costruzione	
Materiale	Stainless steel
Peso	11 111,11 g
Valore nominale	1 mg – 5 kg
Numero	28 pcs

### Compatibile con (Additional Fee)

Pinzette

Pennelli