



More information on the website
mirror.radwag.com/de/info,w1,G3K

Prüfgewichte der M1 Klasse - Knopfgewichte, Set (1 mg - 1 kg), Holz- Etui

OK-401-0083



The drawings, photos and graphics used are for illustrative purposes only.

Technische Daten

| Messtechnische Parameter | |
|--------------------------|-------------|
| OIML-Klasse | M1 |
| Gewicht | 2 111,11 g |
| Nominalwert | 1 mg – 1 kg |
| Anzahl | 25 St. |

Kompatibel mit (Additional Fee)

Pinsel