



More information on the website  
[mirror.radwag.com/en/info,w1,G3K](https://mirror.radwag.com/en/info,w1,G3K)

# M1 Mass Standard, Knob Weights, Set (1 mg - 1 kg), Wooden Box

OK-401-0083



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

### Metrological parameters

OIML Class	M1
Weight	2 111.11 g
Nominal value	1 mg – 1 kg
Quantity	25 pcs

## Compatible with (Additional Fee)

Brushes