



More information on the website
mirror.radwag.com/en/info,w1,GAO

F1 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 g - 1 kg), Plastic Box

OK-202-0058



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Metrological parameters | |
|-------------------------|-----------------|
| OIML Class | F1 |
| Construction | |
| Material | Stainless steel |
| Weight | 2 110 g |
| Nominal value | 1 g – 1 kg |
| Quantity | 13 pcs |

Compatible with (Additional Fee)

Tweezers

Brushes