



More information on the website  
[mirror.radwag.com/de/info,w1,GBN](https://mirror.radwag.com/de/info,w1,GBN)

# Prüfgewicht der F1 Klasse - Zylindergewicht - 2 g

OP-203-0113



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

Messtechnische Parameter	
OIML-Klasse	F1
Konstruktion	
Material	Edelstahl
Nominalwert	2 g

## Kompatibel mit (Additional Fee)

Pinzetten

Pinself