



More information on the website  
[mirror.radwag.com/de/info,w1,GCU](https://mirror.radwag.com/de/info,w1,GCU)

# Prüfgewicht der E1 Klasse - Knopfgewicht - 1 g

OP-503-0046



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

Messtechnische Parameter	
OIML-Klasse	E1
Konstruktion	
Material	Edelstahl
Nominalwert	1 g

## Kompatibel mit (Additional Fee)

Pinzetten

Pinself