



More information on the website  
[mirror.radwag.com/de/info,w1,GLZ](https://mirror.radwag.com/de/info,w1,GLZ)

# Prüfgewicht der E1 Klasse - Knopfgewicht - 50 kg

OP-503-0057



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

Messtechnische Parameter	
OIML-Klasse	E1
Konstruktion	
Material	Edelstahl
Nominalwert	50 kg

## Kompatibel mit (Additional Fee)

Pinsel