



More information on the website
mirror.radwag.com/en/info,w1,GNZ

F2 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 mg - 500 g), Wooden Box



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Metrological parameters | |
|-------------------------|-----------------|
| OIML Class | F2 |
| Construction | |
| Material | Stainless steel |
| Weight | 1 111.11 g |
| Nominal value | 1 mg - 500 g |
| Quantity | 24 pcs |

Compatible with

Tweezers

Brushes