



More information on the website  
[mirror.radwag.com/it/info,w1,GR5](http://mirror.radwag.com/it/info,w1,GR5)

# Pesi standard F2 - cilindrico - 200 g

OP-303-0053



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

### Parametri metrologici

Classe OIML	F2
-------------	----

### Costruzione

Materiale	Stainless steel
-----------	-----------------

Valore nominale	200 g
-----------------	-------

## Compatibile con (Additional Fee)

Pinzette

Pennelli