



More information on the website  
[mirror.radwag.com/us/info,w1,GZV](https://mirror.radwag.com/us/info,w1,GZV)

# SAL/H/PLUS – Anti-vibration table in stainless steel technology

WX-001-0155



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Construction	
Construction	Stainless steel AISI304
Load-bearing capacity	40 kg
Physical parameters	
Packaging dimensions W x D x H	820×711×816 mm
Net weight	61 kg
Gross weight	77 kg
Length	770 mm
Width	680 mm
Height	815±10 mm

## Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration  
 MYA 5Y.P Microbalance for pipettes calibration  
 MYA 5Y.F Microbalance for Filters

MYA 5Y Microbalance  
 PS R1 Precision Balance  
 XA 5Y Analytical Balance

XA 5Y.M.A Microbalance  
PS R2 Precision Balance  
XA 5Y.A Analytical Balance  
AS 5Y Analytical Balance  
PS R2.H Precision Balance  
Analytical Balance for Filters XA 5Y.F

AS R2 PLUS Analytical Balance  
XA 5Y.KO Manual Mass Comparator  
XA 5Y.M Microbalance  
AS R1 PLUS Analytical Balance  
XA 5Y.M.A.S Microbalance for stents weighing

## Device dimensions W x D x H

