



More information on the website  
[mirror.radwag.com/de/info,w1,HAV](https://mirror.radwag.com/de/info,w1,HAV)

# Prüfgewichte der E2 Klasse - Zylindergewichte, Satz (1 g - 50 g), Kunststoff-Etui



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

| Messtechnische Parameter |           |
|--------------------------|-----------|
| OIML-Klasse              | E2        |
| Konstruktion             |           |
| Material                 | Edelstahl |
| Gewicht                  | 110 g     |
| Nominalwert              | 1 – 50 g  |
| Anzahl                   | 8 St.     |