



[radwag.com](http://radwag.com)



## Cable PT0454

WX-013-0240

More information on the website  
[mirror.radwag.com/br/info,w1,HHO](http://mirror.radwag.com/br/info,w1,HHO)



The drawings, photos and graphics used are for illustrative purposes only.

## Visão geral

PT0454 Cable is used to connect the 5Y and CY10 series balance to a communication module Media Box IM02.1. Its standard length is 0,6 m.

## Compatível com (Additional Fee)

Microbalança para calibração de micropipetas XA 5Y.M.A.P  
MediaBox  
Microbalança para calibração de micropipetas MYA 5Y.P  
Balança de precisão 5Y.PM  
Microbalança para filtros MYA 5Y.F  
Microbalança XA 5Y.M.A  
!Terminal de pesagem PUE CY10  
Plataforma Multifuncional CY10  
Balança analítica XA 5Y.A

Ultra-Microbalança UYA 5Y  
Balança de precisão CY10  
Balança analítica para filtros XA 5Y.F  
Microbalança MYA 5Y  
Balança analítica XA 5Y  
Analisador de umidade PMV 5Y  
Microbalança XA 5Y.M  
Ultra-Microbalança para filtros UYA 5Y. F  
Microbalança para pesagem de stents XA 5Y.M.A.S