



More information on the website  
[mirror.radwag.com/br/info,w1,HHO](http://mirror.radwag.com/br/info,w1,HHO)

## Cable PT0454

WX-013-0240



The drawings, photos and graphics used are for illustrative purposes only.

## Visão geral

PT0454 Cable is used to connect the 5Y and CY10 series balance to a communication module Media Box IM02.1. Its standard length is 0,6 m.

## Compatível com (Additional Fee)

Microbalanço para calibrao de micropipetas XA 5Y.M.A.P  
MediaBox  
Microbalanço para calibrao de micropipetas MYA 5Y.P  
Balanço de precisao 5Y.PM  
Microbalanço para filtros MYA 5Y.F  
Microbalanço XA 5Y.M.A  
!Terminal de pesagem PUE CY10  
Plataforma Multifuncional CY10  
Balanço analitica XA 5Y.A

Ultra-Microbalanço UYA 5Y  
Balanço de precisao CY10  
Balanço analitica para filtros XA 5Y.F  
Microbalanço MYA 5Y  
Balanço analitica XA 5Y  
Analisador de umidade PMV 5Y  
Microbalanço XA 5Y.M  
Ultra-Microbalanço para filtros UYA 5Y. F  
Microbalanço para pesagem de stents XA 5Y.M.A.S