



More information on the website
mirror.radwag.com/en/info,w1,HHO

PT0454 Cable

WX-013-0240



The drawings, photos and graphics used are for illustrative purposes only.

Overview

PT0454 Cable is used to connect the 5Y and CY10 series balance to a communication module Media Box IM02.1. Its standard length is 0,6 m.

Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration
MediaBox
MYA 5Y.P Microbalance for pipettes calibration
5Y.PM Precision Balance
MYA 5Y.F Microbalance for Filters
XA 5Y.M.A Microbalance
CY10 Scales for PGC
PUE CY10 Weighing Terminal
CY10 Multifunctional Scale
XA 5Y.A Analytical Balance

UYA 5Y Ultra-Microbalance
CY10 Precision Balance
XA 5Y.F Analytical Balance for Filters
MYA 5Y Microbalance
XA 5Y Analytical Balance
PMV 5Y Moisture Analyzer
XA 5Y.M Microbalance
UYA 5Y.F Ultra-Microbalance for Filters
XA 5Y.M.A.S Microbalance for stents weighing