



More information on the website
mirror.radwag.com/us/info,w1,HHO

PT0454 Cable

WX-013-0240



The drawings, photos and graphics used are for illustrative purposes only.

Overview

RS 232 cable for connecting a balance/scale to a computer.

Cable length: 5 m

Compatible with:

- Scales operated using the PUE C315 terminal
- PUE C315 terminal

Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration
MediaBox
MYA 5Y.P Microbalance for pipettes calibration
5Y.PM Precision balance
MYA 5Y.F Microbalance for Filters
XA 5Y.M.A Microbalance
CY10 Scales for PGC
PUE CY10 Weighing Terminal
CY10 Multifunctional Scale
XA 5Y.A Analytical Balance

UYA 5Y Ultra-Microbalance
CY10 Precision balance
Analytical Balance for Filters XA 5Y.F
MYA 5Y Microbalance
XA 5Y Analytical Balance
PMV 5Y Moisture Analyzer
XA 5Y.M Microbalance
UYA 5Y.F Ultra-Microbalance for Filters
XA 5Y.M.A.S Microbalance for stents weighing