



More information on the website  
[mirror.radwag.com/us/info,w1,HHO](http://mirror.radwag.com/us/info,w1,HHO)

## PT0454 Cable

WX-013-0240



The drawings, photos and graphics used are for illustrative purposes only.

### Overview

RS 232 cable for connecting a balance/scale to a computer.

Cable length: 5 m

Compatible with:

- Scales operated using the PUE C315 terminal
- PUE C315 terminal

### Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration  
MediaBox  
MYA 5Y.P Microbalance for pipettes calibration  
5Y.PM Precision balance  
MYA 5Y.F Microbalance for Filters  
XA 5Y.M.A Microbalance  
CY10 Scales for PGC  
PUE CY10 Weighing Terminal  
CY10 Multifunctional Scale  
XA 5Y.A Analytical Balance

UYA 5Y Ultra-Microbalance  
CY10 Precision balance  
Analytical Balance for Filters XA 5Y.F  
MYA 5Y Microbalance  
XA 5Y Analytical Balance  
PMV 5Y Moisture Analyzer  
XA 5Y.M Microbalance  
UYA 5Y.F Ultra-Microbalance for Filters  
XA 5Y.M.A.S Microbalance for stents weighing