



More information on the website
mirror.radwag.com/de/info,w1,HNE

Prüfgewichte der M1 Klasse - Knopfgewichte, Set (1 mg - 2 kg), Kunststoff-Etui

OK-402-0076



The drawings, photos and graphics used are for illustrative purposes only.

Technische Daten

Messtechnische Parameter	
OIML-Klasse	M1
Gewicht	6 111,11 g
Nominalwert	1 mg – 2 kg
Anzahl	27 St.

Kompatibel mit (Additional Fee)

Pinsel