



More information on the website
mirror.radwag.com/fr/info,w1,HSS

UL-03-1 (AS)



The drawings, photos and graphics used are for illustrative purposes only.

Paramètres Techniques

Paramètres physiques

| | |
|---------------------|---------------|
| Dimensions de colis | 110×102×56 mm |
| Masse nette | 103 g |
| Masse brute | 145 g |

Compatible avec

Balance analytique AS X7
Balances analytiques AS X2 PLUS
Balance analytique AS 5Y

Balance analytique AS R2 PLUS
Balance analytique AS R1 PLUS