



More information on the website
mirror.radwag.com/it/info,w1,HSS

Supporti da laboratorio UL-03-1 (AS)



The drawings, photos and graphics used are for illustrative purposes only.

Specifiche

Parametri fisici

Dimensioni del pacco	110×102×56 mm
Peso netto	103 g
Peso lordo	145 g

Compatibile con

AS X7 Bilancia analitica
Bilance analitiche AS X2 PLUS
Bilancia analitica AS 5Y

Bilancia analitica AS R2 PLUS
Bilancia analitica AS R1 PLUS