



More information on the website  
[mirror.radwag.com/tr/info,w1,HSS](http://mirror.radwag.com/tr/info,w1,HSS)

## UL-03-1 (AS) Laboratuvar için tutucu

WX-001-0032



The drawings, photos and graphics used are for illustrative purposes only.

### Teknik Veriler

#### Fiziksel parametreler

Ambalaj boyutları	110×102×56 mm
Net ağırlık	103 g
Brüt ağırlık	145 g

### İle uyumlu (Additional Fee)

AS X7 Analitik Teraziler  
AS X2 PLUS Analitik Teraziler  
AS 5Y Analitik Teraziler

AS R2 PLUS Analitik Teraziler  
AS R1 PLUS Analitik Teraziler