



More information on the website
mirror.radwag.com/tr/info,w1,HSS

UL-03-1 (AS) Laboratuvar için tutucu

WX-001-0032



The drawings, photos and graphics used are for illustrative purposes only.

Teknik Veriler

Fiziksel parametreler

Ambalaj boyutları	110×102×56 mm
Net ağırlık	103 g
Brüt ağırlık	145 g

İle uyumlu (Additional Fee)

AS X7 Analitik Teraziler
AS X2 PLUS Analitik Teraziler
AS 5Y Analitik Teraziler

AS R2 PLUS Analitik Teraziler
AS R1 PLUS Analitik Teraziler