



More information on the website  
[mirror.radwag.com/de/info,w1,HXW](https://mirror.radwag.com/de/info,w1,HXW)

# Prüfgewicht der E1 Klasse - Plättchenform - 20 mg

OP-503-0041



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

Messtechnische Parameter	
OIML-Klasse	E1
Konstruktion	
Material	Edelstahl
Nominalwert	20 mg

## Kompatibel mit (Additional Fee)

Pinzetten

Pinself