



More information on the website  
[mirror.radwag.com/en/info,w1,l26](https://mirror.radwag.com/en/info,w1,l26)

# Stainless Steel Wire, Set (1 mg - 500 mg), Wooden Box

OK-301-0025



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F2
Construction	
Material	Stainless steel
Weight	1.11 g
Nominal value	1 mg - 500 mg
Weights quantity	12

## Compatible with (Additional Fee)

Tweezers

Brushes