



More information on the website
mirror.radwag.com/it/info,w1,IC2

Cilindrico, set (1 mg - 50 g), cassa di legno

OK-301-0037



The drawings, photos and graphics used are for illustrative purposes only.

Specifiche

| Parametri metrologici | |
|-----------------------|-----------------|
| Classe OIML | F2 |
| Costruzione | |
| Materiale | Stainless steel |
| Peso | 111,11 g |
| Valore nominale | 1 mg – 50 g |
| Numero | 20 pcs |

Compatibile con (Additional Fee)

Pinzette

Pennelli