



More information on the website  
[mirror.radwag.com/en/info,w1,IQ5](http://mirror.radwag.com/en/info,w1,IQ5)

# F2 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 g - 2 kg), Wooden Box

OK-301-0032



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	F2
Construction	
Material	Stainless steel
Weight	6 110 g
Nominal value	1 g – 2 kg
Quantity	15 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes