



More information on the website  
[mirror.radwag.com/it/info,w1,lzv](https://mirror.radwag.com/it/info,w1,lzv)

# Cilindrico, set (1 g - 2 kg), imballaggi in plastica

OK-502-0009



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

Parametri metrologici	
Classe OIML	E1
Costruzione	
Materiale	Stainless steel
Peso	6 110 g
Valore nominale	1 g – 2 kg
Numero	15 pcs

## Compatibile con (Additional Fee)

Pinzette

Pennelli