



More information on the website  
[mirror.radwag.com/en/info,w1,J0Y](https://mirror.radwag.com/en/info,w1,J0Y)

# Knob weights, Set (1 mg - 10 kg), Wooden Box

OK-501-0043



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Metrological parameters	
OIML Class	E1
Construction	
Material	Stainless steel
Weight	21 111.11 g
Nominal value	1 mg – 10 kg
Quantity	29 pcs

## Compatible with (Additional Fee)

Tweezers

Brushes