



More information on the website  
[mirror.radwag.com/de/info,w1,JL7](https://mirror.radwag.com/de/info,w1,JL7)

# Prüfgewicht der E1 Klasse - Knopfgewicht - 20 g

OP-503-0050



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

Messtechnische Parameter	
OIML-Klasse	E1
Konstruktion	
Material	Edelstahl
Nominalwert	20 g

## Kompatibel mit (Additional Fee)

Pinzetten

Pinself