



More information on the website  
[mirror.radwag.com/de/info,w1,JQP](http://mirror.radwag.com/de/info,w1,JQP)

# Prüfgewichte der E1 Klasse - Knopfgewichte, Satz (1 mg - 1 kg), Kunststoff-Etui

OK-502-0018



The drawings, photos and graphics used are for illustrative purposes only.

## Technische Daten

Messtechnische Parameter	
OIML-Klasse	E1
Konstruktion	
Material	Edelstahl
Gewicht	2 111,11 g
Nominalwert	1 mg – 1 kg
Anzahl	25 St.

## Kompatibel mit (Additional Fee)

Pinzetten

Pinzel